

Measurement of GTP-Binding and GTPase Activity of Heterotrimeric G α Proteins

Swarup Roy Choudhury, Corey S. Westfall, Dieter Hackenberg, and Sona Pandey

Abstract

Heterotrimeric G-proteins are important signaling intermediates in all eukaryotes. These proteins link signal perception by a cell surface localized receptor to the downstream effectors of a given signaling pathways. The minimal core of the heterotrimeric G-protein complex consists of G α , G β , and G γ subunits, the G protein coupled receptor (GPCR) and the regulator of G-protein signaling (RGS) proteins. Signal transduction by heterotrimeric G-proteins is controlled by the distinct biochemical activities of G α protein, which binds and hydrolyses GTP. Evaluation of the rate of GTP binding, the rate of GTP hydrolysis, and the rate of GTP/GDP exchange on G α protein are required to better understand the mechanistic aspects of heterotrimeric G-protein signaling, which remains significantly limited for the plant G-proteins. Here we describe the optimized methods for measurement of the distinct biochemical activities of the Arabidopsis G α protein.

Key words Heterotrimeric G-protein, G α protein, GTP-binding, GTP hydrolysis, GTPase activity, GDP/GTP exchange, Arabidopsis GPA1

1 Introduction

G-proteins, localized at the inner surface of the plasma membrane, are key signaling intermediates in eukaryotes [1, 2]. In mammals, the importance of G-proteins in regulating fundamental signaling pathways involved in sensory perception (vision, olfaction, and taste), neurotransmission, hormone perception, and immunity-related cues has prompted in-depth characterization [1–3]. Such studies have revealed an elegant signaling mechanism, where the G α subunit acts as a bimodal molecular switch, alternating between signal-dependent G α ·GDP and G α ·GTP conformations. In the inactive state, G α ·GDP associates with the G $\beta\gamma$ dimer, which represents the “off” signaling status. Signal perception by a G-protein

Swarup Roy Choudhury, Corey S. Westfall and Dieter Hackenberg have contributed equally to this work.

coupled receptor (GPCR) causes the exchange of GTP for GDP on $G\alpha$. The GPCR thus acts as a guanine nucleotide exchange factor (GEF) causing the dissociation of the $G\alpha\beta\gamma$ trimer. The active $G\alpha$ -GTP and the freed $G\beta\gamma$ dimer then transduce the signal downstream by interaction with different effectors. The intrinsic GTPase activity of $G\alpha$ causes hydrolysis of the bound GTP, regenerating $G\alpha$ -GDP, which re-associates with the $G\beta\gamma$ dimer [1, 2]. The GTPase activity of the $G\alpha$ proteins is regulated by a group of accessory proteins, e.g., the regulators of G-protein signaling (RGS) proteins. The RGS proteins, by enhancing the GTPase activity of $G\alpha$, accelerate the rate of G-protein cycle [4, 5]. This classic signaling mechanism thus entails three biochemically distinct reactions of the $G\alpha$ of the heterotrimer: *the rate of GTP binding, the rate of GTP hydrolysis, and the rate of GTP/GDP exchange*, which control all physiological responses regulated by G-proteins.

The biochemical activities of mammalian $G\alpha$ proteins have been characterized in great detail, and it has been established that the rate of GTP-binding on $G\alpha$ proteins is the rate limiting step for mammalian G-protein signaling. In contrast, only few plant $G\alpha$ proteins have been characterized at the biochemical level, and many of them exhibit relatively higher rates of GTP-binding but slower GTPase activities [6–11]. This suggests that the GTPase activity of $G\alpha$ could be the rate limiting step of the G-protein cycle in higher plants, although detailed characterization of multiple $G\alpha$ proteins from different plant lineages is needed to make a generalized statement. In the following sections, we describe optimized methods for purification of plant $G\alpha$ proteins and their biochemical characterization using radiolabeled or fluorescently labeled tags.

2 Materials

2.1 Reagents

1. Purified $G\alpha$ protein (at least 95 % pure).
2. Purified recombinant RGS domain (for performing GTPase assay in the presence of RGS protein).
3. 5 mM GTP (Guanosine triphosphate).
4. 5 mM GDP (Guanosine diphosphate).
5. 5 mM ATP (Adenosine triphosphate).
6. 5 mM ADP (Adenosine diphosphate).
7. [^{35}S]GTP γ S (Perkin-Elmer, 1,250 Ci/mmol).
8. [^{35}S]GTP γ S-binding reaction buffer (Tris pH 8.0, 50 mM; MgCl_2 , 10 mM; DTT, 1 mM).
9. [^{35}S]GTP γ S-binding wash buffer (Tris pH 8.0, 20 mM; MgCl_2 , 25 mM; NaCl, 100 mM).
10. GF/B filter disks (Millipore).
11. Scintillation cocktail (OptiPhase HiSafe 2, Perkin-Elmer).

12. [$\alpha^{32}\text{P}$]GTP (Perkin-Elmer, 3,000 Ci/mmol).
13. PEI cellulose TLC plates (Sigma-Aldrich).
14. TLC assay reaction buffer (Tris pH 8.0, 50 mM; MgCl_2 , 10 mM; EDTA, 1 mM; DTT, 1 mM).
15. TLC assay stop solution (EDTA, pH 8.0, 0.5 M).
16. TLC assay developer solution (KH_2PO_4 , pH 3.4, 0.5 M).
17. BODIPY-assay reaction buffer (Tris pH 8.0, 10 mM; MgCl_2 , 10 mM).
18. 5 mM BODIPY-GTP γ S (Invitrogen).
19. 5 mM BODIPY-GTP FL (Invitrogen).
20. Bovine serum albumin (fraction V powder, Sigma-Aldrich).
21. 5 mM Mant-GTP (Invitrogen).
22. 5 mM Mant-GDP (Invitrogen).
23. Mant reaction Buffer (Tris pH 8.0, 20 mM; NaCl, 100 mM; MgCl_2 , 10 mM).

2.2 Equipment

1. Water bath (30 °C).
2. Vacuum-filtration device.
3. Scintillation vials.
4. Scintillation counter.
5. Cassette and phosphor screen.
6. Phosphorimager.
7. 96 Well Polystyrene Microplates (Greiner Bio One).
8. Fluorescence microplate reader (FLUOstar Optima, BMG Lab Technologies or equivalent).
9. Amicon® Centrifugal Filters (10 KDa cut off).
10. Olis® DM45 spectrofluorimeter with a 150-W xenon lamp, stopped-flow accessory, and water bath.
11. 420 nm cutoff-filter.
12. Fitting software such as KinTek Global Kinetic Explorer [12].

3 Methods

3.1 Radiolabeled GTP-Binding Assay

1. Adjust purified, recombinant $\text{G}\alpha$ protein concentration to 100 nM in 200 μL of [^{35}S]GTP γ S-binding reaction buffer.
2. Add 0.2 μM [^{35}S]GTP γ S and incubate samples at 30 °C in a water bath for up to 2 h (*see Note 1*).
3. Take a small reaction aliquot (5–10 μL) in a separate tube kept chilled on ice at desired time points (*see Note 2*).
4. Add 1 mL of ice-cold wash buffer to the aliquot immediately to stop the reaction. Keep tubes on ice.

5. Wash a GF/B filter with ice cold wash buffer using a vacuum filtration device.
6. Filter the reaction mixture through the GF/B filter, followed by washing with 3–5 mL of ice cold wash buffer.
7. Dry GF/B filters (air dry or oven dry) and place them in a scintillation vial.
8. Add 5 mL of scintillation cocktail to each vial and measure the incorporated radioactivity using a scintillation counter.
9. Analyze and fit data using appropriate software of your choice.

3.2 GTPase Assay Using Thin Layer Chromatography

1. Adjust purified, recombinant $G\alpha$ protein concentration to 100 nM in 200 μ L of TLC assay reaction buffer.
2. Add 20 pmol [$\alpha^{32}P$]GTP and incubate samples at 30 °C in a water bath for up to 2 h.
3. Take a 10 μ L aliquot in a separate tube kept chilled on ice at desired time points (*see* **Notes 1** and **3**).
4. Add 10 μ L of 0.5 M EDTA to the tubes immediately to stop the reaction. Keep tubes on ice.
5. Spot 1 μ L of reaction on to PEI-cellulose TLC plates and air dry. Repeat this step two more times, each time spotting on top of the original spot.
6. Develop plate in 0.5 M KH_2PO_4 (pH 3.4) solution and dry.
7. Expose plate to a phosphorimager screen for 6–12 h. Scan using a phosphorimager.

3.3 GTP Binding Assay Using BODIPY- GTP or BODIPY GTP γ S

1. Use 50 mL BODIPY assay buffer to make 10 nM BODIPY-GTP γ S-FL or BODIPY-GTP-FL solutions. BODIPY-GTP γ S-FL is used for detection of GTP-binding only while BODIPY-GTP-FL can be used for estimation of both GTP-binding and GTP-hydrolysis.
2. Adjust $G\alpha$ protein concentration to 200 nM in BODIPY assay buffer.
3. Add 100 μ L of 200 nM $G\alpha$ protein in 100 μ L assay buffer in 96 well, flat bottom polystyrene microplates and mix properly. A total three to five replicates for each reaction are highly recommended. Reaction mixes containing 100 nM BSA and without any protein should be included as negative controls with each run.
4. Wash the injection tubes getting in and out of the fluorescent plate reader with sterile water, followed by BODIPY-assay buffer.
5. Add 100 μ L of 10 nM BODIPY-GTP γ S or 10 nM BODIPY-GTP FL solutions to the protein solution to initiate the reaction using the automatic injector of the plate reader. The final protein concentration is 100 nM.

6. Record fluorescent reading using following set up for FLUOstar Optima, BMG Lab Technologies fluorescent plate reader: Positioning delay-0.1 s, number of kinetic windows-1, number of cycles-100, number of flashes per cycle-10, filter and integration-fluorescence intensity, number of multichromatics-1, excitation filter-485 nm, emission filter-520, gain-1,500, pump speed-420 $\mu\text{L/s}$, pump prime volume-500 μL , temperature-25 $^{\circ}\text{C}$, data recording-100 cycles every 16 s (total time per read: ~ 30 min). These parameters can be adjusted to desired specifications for any available fluorescent plate reader.
7. For competition experiments, 5 μM of non-labeled nucleotides (GTP/GDP/ATP/ADP) can be added to the assay buffer before starting the reaction. The additives should also be included in the negative controls.
8. For the effect of RGS protein on the rate of GTP hydrolysis, recombinant RGS domain can be added to the reaction mix at a stoichiometric ratio of RGS to $\text{G}\alpha$ at 2:1. A reaction mix containing only RGS protein should be included as a negative control for these assays.
9. Subtract fluorescent readings for the negative controls at each time point and plot relative fluorescence of the proteins as a factor of time using any software of your choice. When using BODIPY GTP-FL, the increase in fluorescent with time denotes GTP-binding and subsequent decrease in fluorescent with time denotes GTP hydrolysis by $\text{G}\alpha$ protein.

3.4 GTP-Binding Using Mant-GTP

1. Using the Amicon $\text{\textcircled{C}}$, exchange purified protein into Mant-GTP reaction buffer and concentrate to 2 μM , keep on ice (*see Note 4*).
2. 5 mM Mant-GTP should be diluted into reaction buffer to form 20, 30, 40, 50, 100, 160 μM , keep on ice (*see Note 5*).
3. Wash fluorometer injection tubes (Olis $\text{\textcircled{C}}$ DM45 spectrofluorimeter or equivalent) 3 times with reaction buffer.
4. Fill stopped-flow syringes with $\text{G}\alpha$ protein solution and Mant-GTP solution.
5. Run 2 injections to wash machine.
6. Let solutions come up to 20 $^{\circ}\text{C}$ for 10 min.
7. Set excitation wavelength to 280 nm and use the >420 nm cut-off for emission. This actually tracks the FRET between a tryptophan and the Mant group. If there is no FRET, an excitation of 350 nm can be used to directly measure Mant fluorescence. FRET is usually less noisy, so it is recommended if possible.
8. Run 1 injection, adjust PMT voltage as needed.
9. Run 9–11 injections, collecting around 1–2 s per injection.
10. Repeat **steps 3–7** for each Mant-GTP concentration.

3.5 Analysis of Mant-GTP Binding Data

1. Normalize all sample runs.
2. Average normalized samples.
3. Repeat for each Mant-GTP concentrations.
4. Depending on the protein, fit the data to either a single or double exponential curve.
5. A graph of Mant-GTP versus observed rate constants should be a line with a slope equal to the binding constant and an x -intercept of the dissociation constant. The Mant-GTP concentration used should be half of starting concentration due to the 1:1 mixing with the protein solution in the stopped-flow.
6. For data that is fit to double exponential curves, one set of observed rate constants will behave linearly with a slope equal to the binding constant and an x -intercept of the dissociation constant. The second observed rate constant should be constant for all Mant-GTP concentrations and corresponds to a conformational change upon binding [13] (*see Note 6*). The equation for the linear fit is $k_{\text{obs}} = k_{-1} + k_1[\text{Ligand}]$.

3.6 Measurement of GDP/GTP Exchange

1. Using the Amicon®, exchange purified protein into reaction buffer and concentrate to 2 μM , keep on ice as in Subheading 3.4 step 1 (*see Note 4*).
2. Add Mant-GDP directly to the protein solution to a final concentration of 5 μM (*see Note 7*).
3. Make solutions of GTP (no Mant group) of 40, 80, 160, 320 μM in Mant reaction buffer.
4. Fill stopped-flow syringes with $G\alpha$ protein plus Mant-GDP solution and GTP solution.
5. Run 2 injections to wash machine.
6. Let solutions come up to 20 °C for 10 min.
7. Set excitation wavelength to 280 nm and use the >420 nm cutoff for emission.
8. Run 1 injection, adjust PMT voltage.
9. Run 9–11 injections collecting 30–45 s per injection.
10. Repeat steps 4–9 for each GTP concentration.

3.7 Analysis of Nucleotide Exchange Data

1. Normalize all 9–11 sample runs and average normalized samples.
2. Repeat for each GTP concentrations.
3. Data for each GTP concentration should be able to fit to a single-exponential curve. If not, removal of beginning points or end points might be necessary.
4. A graph of the observed rate constants should be hyperbolic. If the rate constants do not reach a plateau, higher GTP concentrations might need to be used.

5. The observed rate constant for saturated GTP corresponds to the dissociation rate of the Mant-GDP [14] (*see* **Note 8**).

4 Notes

1. Strictly follow all safety procedures for handling and disposal of radioisotopes and liquid and solid waste containing radiochemicals. Survey the usage area and equipment before and after the use of radioisotopes. Follow your institutional guidelines for recording of radioisotope usage and disposal.
2. The aliquots should be collected every 30 s to 1 min, if possible to achieve better resolution of binding kinetics at the initial time points. Three to five aliquots should be collected for each time point. For strong GTP-binding proteins, 30–60 min incubation is sufficient.
3. The aliquots should be collected every 30 s to 1 min at the initial time points and should be added directly to ice cold tubes containing chilled EDTA. Three to five aliquots should be collected for each time point. Maximum 18 samples can be spotted on a standard TLC plate for each run.
4. Protein concentration needs to be calculated using accurate techniques, such as absorbance at 280 nm.
5. Exact concentrations of protein, Mant-GTP, and Mant-GDP will depend on the protein being used. To aid in analysis, always maintain pseudo-first order kinetics by keeping the nucleotide 10× higher than the protein.
6. This analysis only works if you always maintain pseudo-first order kinetics. If you are not pseudo-first order, you must either fit computationally, using a program like KinTek®, or use a more complicated fitting equation [13].
7. Mant-GDP concentration needs to be high enough to saturate binding sites.
8. Some proteins might have a negative hyperbolic curve for the exchange reaction. *See* ref. 14 for this circumstance.

Acknowledgments

The research in corresponding author's laboratory was supported by a US Department of Agriculture/Agriculture and Food Research Initiative grant (2010-65116-20454) and a National Science Foundation grant (MCB-1157944). C.S.W. was supported by a US Department of Agriculture/Agriculture and Food Research Initiative predoctoral research fellowship (MOW-2010-05240).

References

1. Offermanns S (2003) G-proteins as transducers in transmembrane signalling. *Prog Biophys Mol Biol* 83(2):101–130
2. Cabrera-Vera TM et al (2003) Insights into G protein structure, function, and regulation. *Endocr Rev* 24(6):765–781
3. Oldham WM, Hamm HE (2008) Heterotrimeric G protein activation by G-protein-coupled receptors. *Nat Rev Mol Cell Biol* 9(1):60–71
4. Dohlman HG, Thorner J (1997) RGS proteins and signaling by heterotrimeric G proteins. *J Biol Chem* 272(7):3871–3874
5. Siderovski DP, Willard FS (2005) The GAPs, GEFs, and GDIs of heterotrimeric G-protein alpha subunits. *Int J Biol Sci* 1(2): 51–66
6. Jones JC, Temple BR, Jones AM, Dohlman HG (2011) Functional reconstitution of an atypical G protein heterotrimer and regulator of G protein signaling protein (RGS1) from *Arabidopsis thaliana*. *J Biol Chem* 286(15): 13143–13150
7. Jones JC et al (2011) The crystal structure of a self-activating G protein alpha subunit reveals its distinct mechanism of signal initiation. *Sci Signal* 4(159):ra8
8. Johnston CA et al (2007) GTPase acceleration as the rate-limiting step in Arabidopsis G protein-coupled sugar signaling. *Proc Natl Acad Sci USA* 104(44):17317–17322
9. Seo HS, Choi CH, Lee SY, Cho MJ, Bahk JD (1997) Biochemical characteristics of a rice (*Oryza sativa* L., IR36) G-protein alpha-subunit expressed in *Escherichia coli*. *Biochem J* 324(Pt 1):273–281
10. Bisht NC, Jez JM, & Pandey S (2011) An elaborate heterotrimeric G-protein family from soybean expands the diversity of plant G-protein networks. *New phytol* 190:35–48
11. Roy Choudhury S et al (2012) Two chimeric regulators of G-protein signaling (RGS) proteins differentially modulate soybean heterotrimeric G-protein cycle. *J Biol Chem* 287(21): 17870–17881
12. Johnson KA (2009) Fitting enzyme kinetic data with KinTek Global Kinetic Explorer. *Methods Enzymol* 467:601–626
13. Kozlov AG, Lohman TM (2002) Stopped-flow studies of the kinetics of single-stranded DNA binding and wrapping around the *Escherichia coli* SSB tetramer. *Biochemistry* 41(19): 6032–6044
14. Wu XM, Gutfreund H, Chock PB (1992) Kinetic method for differentiating mechanisms for ligand exchange reactions: application to test for substrate channeling in glycolysis. *Biochemistry* 31(7):2123–2128